



**Diamond**



**Carbide**



**Polishing**



**Sets**

**Dentist**

























***Dumont***

**INSTRUMENTS & CO**

**Quality since 1978**

# Icons

- |  |   |
|--|---|
|  Cavity preparation            |  Angle                           |
|  Crown preparation           |  Rounded edges                 |
|  Removal of old fillings     |  Non cutting tip               |
|  Working on fillings         |  End cutting only              |
|  Root canal preparation      |  Further information available |
|  Oral surgery                |  Non cutting tip               |
|  Crown and bridge technique  |  Crown removal                 |
|  Milling technique           |  Orthodontics                  |
|  Model casting technique     |  Diamond interspersed          |
|  Adjusting occlusal surfaces |  Prophylaxis                   |
|  Rounded tip                 |  Implantology                  |

# READING KEYS

## Diamond coated burs

1. Grit  
2. Shape  
3. Head  $\varnothing$   
4. Shank

Example of an order

3 846 018 3

1 2 3 4

### Grit

● 80 $\mu$ (ISO 524)	MEDIUM
● 110 $\mu$ (ISO 534)	COARSE
● 165 $\mu$ (ISO 544)	VERY COARSE

### Grit

● 60 $\mu$ (ISO 514)	FINE
● 20 $\mu$ (ISO 504)	VERY FINE
● 10 $\mu$ (ISO 494)	ULTRA FINE

HP SHANK Code 1 =  $\varnothing$  2,35 MM



RA SHANK Code 2 =  $\varnothing$  2,35 MM



FG SHANK Code 3 =  $\varnothing$  1,60 MM



## Carbide instruments

1. Shape  
2. Head  $\varnothing$   
3. Shank

Example of an order

C1 - 008 - 1

1 2 3

HP SHANK Code 1 =  $\varnothing$  2,35 MM



RA SHANK Code 2 =  $\varnothing$  2,35 MM



FG SHANK Code 3 =  $\varnothing$  1,60 MM



## Polishing instruments

HP SHANK Code 1 =  $\varnothing$  2,35 MM



RA SHANK Code 2 =  $\varnothing$  2,35 MM



FG SHANK Code 3 =  $\varnothing$  1,60 MM



# Diamond instruments





# Diamond coated burs


## DIAMOND



### 368

**Bud pointed**

806 314 257



Shank

L mm 3,5 3,5 3,5 4,5 4,5 5,0



Order N° Shank

FG ●	524 ...	368-	010		016	018	021	023	-3
FG ●	544 ...	5368-		012		018	021	023	-3
FG ●	534 ...	3368-	010	012	016	018	021	023	-3
FG ●	514 ...	8368-	010		016	018	021	023	-3
FG ●	504 ...	4368-			016	018	021	023	-3
FG ●	494 ...	368UF-						023	-3



### 370

**Occlusal reduction**

806 314 507



Shank

L mm 5,7



Order N° Shank


FG ●	534 ...	3370-	023	-3
------	---------	-------	-----	----



### 379

**Egg/Football**

806 314 263



Shank

L mm 3,0 3,0 3,5 3,5 4,0 4,5



Order N° Shank

FG ●	524 ...	379-	012	014	016	018		023	-3
FG ●	544 ...	5379-		014	016	018	021	023	-3
FG ●	534 ...	3379-	012	014	016	018		023	-3
FG ●	514 ...	8379-	012	014	016	018	021	023	-3
FG ●	504 ...	4379-	012	014	016	018	021	023	-3


Match C44E series



# 390

## Grenade

806 314 274



Shank

L mm 3,5



Order N° Shank

FG ●	524 ...	390-	016	-3
FG ●	514 ...	8390-	016	-3
FG ●	504 ...	4390-	016	-3



# 393

## Interproximal

806 314 465



Shank

L mm 5,0

Order N° Shank


FG ●	524 ...	393-	016	-3
FG ●	534 ...	3393-	016	-3
FG ●	514 ...	8393-	016	-3



# 801

## Round

806 314 001



Shank



Order N° Shank

FG ●	524 ...	801-	006	008	009	010	012	014	016	018	021	023	025	027	029	033	035	-3
FG ●	544 ...	5801-		008	009	010	012	014	016	018	021	023				033	035	-3
FG ●	534 ...	3801-		008	009	010	012	014	016	018	021	023	025		029	033		-3
FG ●	514 ...	8801-		008		010	012	014	016	018	021	023	025		029	033		-3
FG ●	504 ...	4801-		008	009	010	012	014	016	018	021	023			029			-3

806 204 001


RA ●	524 ...	801-		008	009	010	012	014	016	018		023		027				-2
RA ●	534 ...	3801-					012	014										-2



## 801L

### Long Round

806 314 697



Shank




		Order N°					Shank
FG ●	524 ...	801L-	012	014	016	-3	
FG ●	534 ...	3801L-	012	014	016	-3	
FG ●	514 ...	8801L-			016	-3	



## 802

### Round with collar

806 314 002



Shank



		Order N°										Shank
FG ●	524 ...	802-	008	010	012	014	016	018	021	023	-3	
FG ●	544 ...	5802-			012	014	016	018	021	023	-3	
FG ●	534 ...	3802-		010	012	014	016	018	021	023	-3	
FG ●	514 ...	8802-		010	012	014	016	018	021	023	-3	
FG ●	504 ...	4802-	008	010	012	014	016	018	021		-3	



## 802L

### Long Round with collar

806 314 494



Shank




		Order N°			Shank
FG ●	524 ...	802L-	16	-3	
FG ●	534 ...	3802L-	16	-3	
FG ●	514 ...	8802L-	16	-3	



## 803

### Round with long collar

806 314 003



Shank



		Order N°							Shank
FG ●	524 ...	803-	014	016	018	021	023	-3	
FG ●	544 ...	5803-		016	018	021	023	-3	
FG ●	534 ...	3803-	014	016	018	021		-3	
FG ●	514 ...	8803-		016	018	021		-3	
FG ●	504 ...	4803-		016		021		-3	






## 805

### Inverted cone

806 314 010



Shank

L mm	1,0	1,0	1,2	1,4	1,4	1,5	1,9
------	-----	-----	-----	-----	-----	-----	-----



Order N°	Shank								
----------	-------	--	--	--	--	--	--	--	--

FG ●	524 ...	805-	008	010	012	014	016	018		-3
FG ●	544 ...	5805-		010	012	014	016	018	021	-3
FG ●	534 ...	3805-	008	010	012	014	016	018	021	-3
FG ●	514 ...	8805-	008	010	012	014	016	018		-3
FG ●	504 ...	4805-	008	010	012	014	016			-3

806 204 010


RA ●	524...	805-	008	010	012	014	016	018	021	-2
------	--------	------	-----	-----	-----	-----	-----	-----	-----	----



## 806

### Diabolo

806 314 019



Shank

L mm	2,5	2,5	2,5	2,5	2,5	2,5
------	-----	-----	-----	-----	-----	-----



Order N°	Shank								
----------	-------	--	--	--	--	--	--	--	--


FG ●	524 ...	806-	008	010	012	014	016	018		-3
FG ●	544 ...	5806-			012	014	016	018		-3
FG ●	534 ...	3806-	008	010	012	014	016	018		-3
FG ●	514 ...	8806-		010	012	014	016	018		-3
FG ●	504 ...	4806-		010	012	014	016			-3



## 807

### Long inverted cone

806 314 225



Shank

L mm	3,0	3,5	3,5	4,0	5,0	5,0	6,0
------	-----	-----	-----	-----	-----	-----	-----



Order N°	Shank								
----------	-------	--	--	--	--	--	--	--	--

FG ●	524 ...	807-	010	012	014	016	018	021	023	-3
FG ●	544 ...	5807-				016			023	-3
FG ●	534 ...	3807-		012	014	016	018	021	023	-3
FG ●	514 ...	8807-					018		023	-3
FG ●	504 ...	4807-				016				-3



## 807L

### Longer Inverted cone

806 314 226




Shank		L mm		
		7,0		
		Order N°		Shank
FG ●	524...	807L-	018	-3



## 811

### Occlusal/palatal reduction

806 314 038




Shank		L mm				
		4,2	7,0	7,0		
		Order N°			Shank	
FG ●	524 ...	811-	033	037	047	-3
FG ●	544 ...	5811-	033	037	047	-3
FG ●	534 ...	3811-	033	037	037	-3



## 813

### Diabolo

806 314 035




Shank		L mm						
		1,5	1,5	1,5	2,0	2,0		
		Order N°					Shank	
FG ●	524 ...	813-	010	012	014	016	018	-3
FG ●	544 ...	5813-		012	014	016	018	-3
FG ●	534 ...	3813-	010	012	014	016	018	-3
FG ●	514 ...	8813-	010	012	014	016	018	-3
FG ●	504 ...	4813-	010	012	014	016		-3



## 815

### Wheel

806 314 040



Shank		L mm			
		0,7	0,8		
		Order N°		Shank	
FG ●	524 ...	815-	016	025	-3
FG ●	544 ...	5815-		025	-3
FG ●	534 ...	3815-	016	025	-3

# 818

## Wheel

806 314 042



### Shank

L mm	0,7
------	-----




Order N°	Shank
----------	-------

FG ●	524 ...	818-	035	-3
FG ●	544 ...	5818-	035	-3
FG ●	534 ...	3818-	035	-3

# 830

## Pear

806 314 233



### Shank

L mm	2,7	2,7	2,7	2,7	2,7	2,7	5,1
------	-----	-----	-----	-----	-----	-----	-----




Order N°	Shank
----------	-------

FG ●	524 ...	830-	009	010	012	014	016	018		-3
FG ●	544 ...	5830-			012	014	016	018	021	-3
FG ●	534 ...	3830-		010	012	014	016	018		-3
FG ●	514 ...	8830-		010	012	014	016	018		-3
FG ●	504 ...	4830-		010	012	014	016			-3

# 830L

## Long pear

806 314 234



### Shank

L mm	4,0	4,0	5,0	5,0	5,0	5,0
------	-----	-----	-----	-----	-----	-----



Order N°	Shank
----------	-------

FG ●	524 ...	830L-	010	012	014	016	018		-3
FG ●	544 ...	5830L-			014	016	018		-3
FG ●	534 ...	3830L-		012	014	016	018	021	-3
FG ●	514 ...	8830L-		012	014	016	018		-3
FG ●	504 ...	4830L-		012	014				-3

# 833

## Occlusal/Cavity margin trimming

806 314 466



### Shank

L mm	3,6
------	-----



Order N°	Shank
----------	-------


FG ●	524 ...	833-	031	-3
FG ●	544 ...	5833-	031	-3
FG ●	534 ...	3833-	031	-3
FG ●	514 ...	8833-	031	-3



## 834

Depth marker, veneer technique

806 314 552




Shank		L mm	6,0	6,0	6,0
		T mm	0,3	0,4	0,5
		Order N°	Shank		
FG ●	524 ...	834	016	018	021 -3



## 835

Cylinder flat end

806 314 107



Shank		L mm	3,0	4,0	4,0	4,0	4,0	4,0	4,0
		Order N°	Shank						
FG ●	524 ...	835-	008	010	012	014	016	018	-3
FG ●	544 ...	5835-		010	012	014	016	018	-3
FG ●	534 ...	3835-	008	010	012	014	016	018	-3
FG ●	514 ...	8835-		010	012	014			-3
FG ●	504 ...	4835-	008	010	012	014	016		-3

806 204 107


RA ●	524...	835-	008	010	012		016		-2
------	--------	------	-----	-----	-----	--	-----	--	----



## 835KR

Cylinder, rounded edge

806 314 156




Shank		L mm	4,0	4,0	4,0	4,0
		Order N°	Shank			
FG ●	524 ...	835KR-	010	012	014	-3
FG ●	544 ...	5835KR-			014	-3
FG ●	534 ...	3835KR-	010	012	014	016 -3
FG ●	514 ...	8835KR-		012		-3

# 836

## Short parallel shoulder

806 314 110



Shank

L mm	6,0	6,0	6,0	6,0	6,0	6,0
------	-----	-----	-----	-----	-----	-----



Order N°	Shank							
----------	-------	--	--	--	--	--	--	--

FG ●	524 ...	836-	010	012	014	016	018		-3
FG ●	544 ...	5836-	010	012	014	016	018	021	-3
FG ●	534 ...	3836-	010	012	014	016	018		-3
FG ●	514 ...	8836-	010	012	014				-3
FG ●	504 ...	4836-	010	012	014	016			-3


806 204 110

RA ●	524...	836-		012					-2
------	--------	------	--	-----	--	--	--	--	----

# 836KR

## Short parallel shoulder, rounded edge

806 314 157



Shank

L mm	6,0	6,0	6,0	6,0
------	-----	-----	-----	-----




Order N°	Shank			
----------	-------	--	--	--

FG ●	524 ...	836KR-		012				-3
FG ●	544 ...	5836KR-			014			-3
FG ●	534 ...	3836KR-	010	012	014	016		-3
FG ●	514 ...	8836KR-	010	012				-3

# 837

## Parallel shoulder

806 314 111



Shank

L mm	8,0	8,0	8,0	8,0	8,0
------	-----	-----	-----	-----	-----



Order N°	Shank				
----------	-------	--	--	--	--

FG	524 ...	837-	010	012	014			-3
FG ●	544 ...	5837-		012	014	016	018	-3
FG ●	534 ...	3837-	010	012	014	016		-3
FG ●	514 ...	8837-		012	014			-3
FG ●	504 ...	4837-	010	012	014	016		-3

806 204 111


RA ●	524...	837-		012				-2
------	--------	------	--	-----	--	--	--	----



## 837KR

Parallel shoulder,  
rounded edge

806 314 158



Shank		L mm					
		8,0	8,0	8,0	8,0	Shank	
		Order N°				Shank	
FG ●	524 ...	837KR-		012	014		-3
FG ●	544 ...	5837KR-		012		016	-3
FG ●	534 ...	3837KR-	010	012	014	016	-3
FG ●	514 ...	8837KR-	010	012	014	016	-3


Match C244KR



## 837L

Long parallel shoulder

806 314 112




Shank		L mm					
		10,0	10,0	10,0	10,0	Shank	
		Order N°				Shank	
FG ●	524 ...	837L-		012			-3
FG ●	534 ...	3837L-	010	012	014	016	-3



## 838

Short cylinder rounded edge

806 314 137




Shank		L mm							
		3,0	4,0	4,0	4,0	4,0	4,0	Shank	
		Order N°						Shank	
FG ●	524 ...	838-		010	012	014	016	018	-3
FG ●	544 ...	5838-			012	014	016	018	-3
FG ●	534 ...	3838-	008	010	012	014	016	018	-3
FG ●	514 ...	8838-	008	010	012	014	016	018	-3
FG ●	504 ...	4838-	008	010	012	014	016		-3



## 839

Margin refinement

806 314 150




Shank		L mm			
		Order N°		Shank	
FG ●	524 ...	839	012	014	-3

For lowering the preparation limit after shoulder preparation with 837/837L or flat preparation floor inside cavity.

# 840

## Margin refinement

806 314 722



Shank



L mm

Order N°

Shank


FG ●	524 ...	840	012	014	016	-3
------	---------	-----	-----	-----	-----	----

For lowering the preparation limit after shoulder preparation with 837KR or flat preparation floor inside cavity.

# 845

## Flat end taper

806 314 168



Shank



L mm

3,0 4,0 4,0 4,0

Order N°

Shank

FG ●	524 ...	845-	008	010		014	-3
FG ●	534 ...	3845-	008	010		014	-3
FG ●	514 ...	8845-		010		014	-3


806 204 168

RA ●	524...	845-	008	010	012		-2
------	--------	------	-----	-----	-----	--	----

# 845KR

## Flat end taper, rounded edge

806 314 544



Shank



L mm

4,0 4,0 4,0 4,0

Order N°

Shank

FG ●	524 ...	845KR-	016	018	021	025	-3
FG ●	534 ...	3845KR-	016	018		025	-3
FG ●	514 ...	8845KR-	016	018		025	-3
FG ●	504 ...	4845KR-	016	018		025	-3






## 846

### Short tapered shoulder

806 314 171



Shank		L mm						
		6,0	6,0	6,0	6,0	6,0		
		Order N°					Shank	
FG●	524 ...	846-	012		016	018	021	-3
FG●	544 ...	5846-				018	021	-3
FG●	534 ...	3846-	012	014	016	018		-3
FG●	514 ...	8846-	012	014	016	018		-3
FG●	504 ...	4846-			016			-3

## 846KR

### Short tapered shoulder, rounded edge

806 314 545




Shank		L mm			
		6,0			
		Order N°		Shank	
FG●	524 ...	846KR-	016	-3	
FG●	514 ...	8846KR-	016	-3	
FG●	504 ...	4846KR-	016	-3	



## 847

### Tapered shoulder

806 314 172




Shank		L mm						
		8,0	8,0	8,0	8,0	8,0		
		Order N°					Shank	
FG●	524 ...	847-	008	012	014	016	018	-3
FG●	544 ...	5847-			014	016	018	-3
FG●	534 ...	3847-		012	014	016		-3
FG●	514 ...	8847-	008	012	014	016		-3
FG●	504 ...	4847-		012	014			-3



## 847KR

### Tapered shoulder, rounded edge

806 314 546



Shank		L mm				
		8,0	8,0	8,0		
		Order N°			Shank	
FG●	524 ...	847KR-	016		023	-3
FG●	534 ...	3847KR-	016	018		-3
FG●	514 ...	8847KR-	016	018	023	-3
FG●	504 ...	4847KR-	016			-3






# 848

## Long tapered shoulder

806 314 173



Shank

L mm 10,0 10,0 10,0 10,0 10,0 10,0



Order N° Shank

FG ●	524 ...	848-	010	012	014	016	018	023	-3
FG ●	544 ...	5848-		012		016	018	023	-3
FG ●	534 ...	3848-	010	012	014	016	018	023	-3
FG ●	514 ...	8848-	010	012		016	018	023	-3
FG ●	504 ...	4848-	010		014	016			-3



# 848KR

## Long tapered shoulder, rounded edge

806 314 553



Shank

L mm 10,0



Order N° Shank


FG ●	524 ...	848KR-	016	-3
FG ●	514 ...	8848KR-	016	-3



# 848L

## Extra long tapered shoulder

806 314 174



Shank

L mm 11,5 11,5 11,5



Order N° Shank


FG ●	534 ...	3848L-	014	016	018	-3
------	---------	--------	-----	-----	-----	----



# 849

## Tapered round

806 314 194



Shank

L mm 4,0 4,0 4,0 4,0



Order N° Shank


FG ●	524 ...	849-	010	014	016	018	-3
FG ●	544 ...	5849-			016	018	-3
FG ●	534 ...	3849-		014	016	018	-3
FG ●	514 ...	8849-		014	016	018	-3
FG ●	504 ...	4849-		014	016		-3



## 850

### Tapered round, long

806 314 199



Shank

L mm 5,0 5,0 5,0 5,0 5,0



Order N° Shank


FG ●	524 ...	850-	012	014	016	018	031	-3
FG ●	544 ...	5850-		014	016	018	031	-3
FG ●	534 ...	3850-		014	016	018	031	-3
FG ●	514 ...	8850-	012	014	016	018		-3
FG ●	504 ...	4850-		014	016			-3



## 851

### Safe end, round end tapered

806 314 219



Shank

L mm 8,0 8,0



Order N° Shank


FG ●	524 ...	851	010	012	-3
------	---------	-----	-----	-----	----



## 852

### Short pointed

806 314 164



Shank

L mm 6,0 6,0 6,0 6,0 6,0 6,0 6,0 6,0 6,0



Order N° Shank

FG ●	524 ...	852-	008	010	012	014	016		023		045	-3
FG ●	544 ...	5852-		010	012	014	016	018	023		045	-3
FG ●	534 ...	3852-	008	010	012	014	016		023		045	-3
FG ●	514 ...	8852-	008	010	012	014	016		023		045	-3
FG ●	504 ...	4852-	008	010	012	014						-3
FG ●	494 ...	852UF-			012	014						-3

806 204 164


RA ●	524 ...	852-								037		-2
------	---------	------	--	--	--	--	--	--	--	-----	--	----

Match C333

# 855

## Short tapered chamfer, round

806 314 197



Shank

L mm 6,0 6,0 6,0 6,0 6,0 7,0




Order N° Shank

FG ●	524 ...	855-	010	012	014	016		025	-3
FG ●	544 ...	5855-			014		018		-3
FG ●	534 ...	3855-		012	014				-3
FG ●	514 ...	8855-		012					-3
FG ●	504 ...	4855-		012					-3

# 856

## Tapered chamfer, round

806 314 198



Shank

L mm 8,0 8,0 8,0 8,0 8,0 8,0 8,0 8,0



Order N° Shank


FG ●	524 ...	856-	010	012	014	016	018	021	023	033	-3
FG ●	544 ...	5856-		012	014	016	018	021	023	033	-3
FG ●	534 ...	3856-	010	012	014	016	018	021	023	033	-3
FG ●	514 ...	8856-	010	012	014	016	018	021	023	033	-3
FG ●	504 ...	4856-	010	012	014	016		021			-3

Match C375LR series

# 856L

## Long tapered chamfer, round

806 314 199



Shank

L mm 9,0 9,0 9,0



Order N° Shank


FG ●	524 ...	856L-	014	016		-3
FG ●	544 ...	5856L-	014	016	018	-3
FG ●	534 ...	3856L-	014	016	018	-3
FG ●	514 ...	8856L-	014	016	018	-3
FG ●	504 ...	4856L-	014			-3



# 857

## Safe end, round end tapered

806 314 220




Shank		L mm 10,0 10,0 10,0					
		Order N°				Shank	
FG ●	524 ...	857-	014	016	018		-3
FG ●	544 ...	5857-		016	018		-3
FG ●	534 ...	3857-	014	016	018		-3
FG ●	514 ...	8857-	014	016	018		-3



# 858

## Pointed

806 314 165




Shank		L mm 8,0 8,0 8,0					
		Order N°				Shank	
FG ●	524 ...	858-	010	012	014		-3
FG ●	544 ...	5858-		012	014		-3
FG ●	534 ...	3858-	010	012	014		-3
FG ●	514 ...	8858-	010	012	014		-3
FG ●	504 ...	4858-		012	014		-3



# 859

## Needle

806 314 166



Shank		L mm 10,0 10,0 10,0 10,0 10,0 10,0								
		Order N°							Shank	
FG ●	524 ...	859-	010	012	014	016	018			-3
FG ●	544 ...	5859-		012	014	016	018	021		-3
FG ●	534 ...	3859-	010	012	014	016	018			-3
FG ●	514 ...	8859-	010	012	014					-3
FG ●	504 ...	4859-	010	012	014	016				-3
FG ●	494 ...	859UF-			014	016				-3

Match C335

# 859L

## Long needle

806 314 167



### Shank



L mm 11,0


Order N° Shank

FG ●	524 ...	859L-	014	-3
FG ●	534 ...	3859L-	014	-3
FG ●	514 ...	8859L-	014	-3

# 860

## Short flame

806 314 245



### Shank



L mm 4,0 5,0 5,0 5,0

Order N° Shank

FG ●	524 ...	860-			014	016	-3
FG ●	544 ...	5860-				016	-3
FG ●	534 ...	3860-	010	012	014	016	-3
FG ●	514 ...	8860-	010	012	014	016	-3
FG ●	504 ...	4860-	010	012	014		-3
FG ●	494 ...	860UF-		012			-3

# 861

## Flame

806 314 248



### Shank



L mm 6,0 6,0

Order N° Shank

FG ●	524 ...	861-	012		-3
FG ●	544 ...	5861-		014	-3
FG ●	534 ...	3861-	012	014	-3
FG ●	514 ...	8861-	012		-3






# 862

## Flame

806 314 249



### Shank

L mm	8,0	8,0	8,0	8,0	8,0
------	-----	-----	-----	-----	-----



Order N°	Shank						
----------	-------	--	--	--	--	--	--

FG ●	524 ...	862-	010	012	014	016		-3
FG ●	544 ...	5862-		012	014	016	018	-3
FG ●	534 ...	3862-	010	012	014	016	018	-3
FG ●	514 ...	8862-	010	012	014	016		-3
FG ●	504 ...	4862-	010	012	014	016		-3


Match C48L



# 863

## Long flame

806 314 250



### Shank

L mm	10,0	10,0	10,0	10,0
------	------	------	------	------



Order N°	Shank			
----------	-------	--	--	--

FG ●	524 ...	863-	012	014	016	018		-3
FG ●	544 ...	5863-	012	014	016	018		-3
FG ●	534 ...	3863-	012	014	016			-3
FG ●	514 ...	8863-	012	014	016			-3
FG ●	504 ...	4863-	012	014				-3



# 863GK

## Long flame, safe end

806 314 256



### Shank

L mm	10,0
------	------



Order N°	Shank
----------	-------

FG ●	524 ...	863GK-	012	-3
------	---------	--------	-----	----



# 864

## Extra long flame

806 314 251



### Shank

L mm	12,0
------	------




Order N°	Shank
----------	-------

FG ●	524 ...	864-	012	-3
FG ●	534 ...	3864-	012	-3
FG ●	514 ...	8864-	012	-3

# 868

## Tapered round

806 314 223



Shank

L mm 8,0 8,0 8,0 8,0



Order N° Shank

FG ●	524 ...	868-	012	014	016		-3
FG ●	544 ...	5868-		014	016	018	-3
FG ●	534 ...	3868-	012	014	016		-3
FG ●	514 ...	8868-	012	014	016		-3
FG ●	504 ...	4868-	012	014			-3

# 874

## Adjusting preparation limit

806 314 534



Shank

L mm 2,0



Order N° Shank

FG ●	524 ...	874-	009	-3
FG ●	534 ...	3874-	009	-3
FG ●	514 ...	8874-	009	-3

# 875

## Adjusting preparation limit

806 314 535



Shank

L mm 3,0



Order N° Shank

FG ●	524 ...	875-	009	-3
FG ●	534 ...	3875-	009	-3
FG ●	514 ...	8875-	009	-3

# 876

## Short parallel chamfer, torpedo

806 314 287



Shank

L mm 5,0



Order N° Shank


FG ●	524 ...	876-	010	-3
FG ●	534 ...	3876-	010	-3
FG ●	514 ...	8876-	010	-3



## 877

### Short parallel chamfer, torpedo

806 314 288



Shank		L mm					
		6,0	6,0	6,0	6,0		
		Order N°				Shank	
FG ●	524 ...	877-	010	012	014	016	-3
FG ●	544 ...	5877-		012			-3
FG ●	534 ...	3877-	010	012		016	-3
FG ●	514 ...	8877-	010	012		016	-3
FG ●	504 ...	4877-	010				-3


Match C243



## 877K

### Short tapered chamfer, torpedo

806 314 297



Shank		L mm			
		6,0	6,0		
		Order N°		Shank	
FG ●	524 ...	877K-		016	-3
FG ●	544 ...	5877K-		016	-3
FG ●	534 ...	3877K-	014	016	-3
FG ●	514 ...	8877K-	014	016	-3


Match C282K



## 878

### Parallel chamfer, torpedo

806 314 289



Shank		L mm					
		8,0	8,0	8,0	8,0		
		Order N°				Shank	
FG ●	524 ...	878-	010	012	014	016	-3
FG ●	544 ...	5878-		012			-3
FG ●	534 ...	3878-	010	012		016	-3
FG ●	514 ...	8878-	010	012		016	-3
FG ●	504 ...	4878-	010				-3

Match C283






## 878K

### Tapered chamfer, torpedo

806 314 298



Shank		L mm								
		8,0	8,0	8,0	8,0	8,0	8,0	Shank		
		Order N°							Shank	
FG ●	524 ...	878K-	012	014	016	018	021		-3	
FG ●	544 ...	5878K-	012	014	016	018	021	023	-3	
FG ●	534 ...	3878K-	012	014	016	018	021		-3	
FG ●	514 ...	8878K-	012	014	016	018	021		-3	
FG ●	504 ...	4878K-	012	014	016				-3	


Match C244K



## 879

### Long parallel chamfer, torpedo

806 314 290




Shank		L mm								
		10,0	10,0	10,0	10,0	Shank				
		Order N°					Shank			
FG ●	524 ...	879-	012	014	016					-3
FG ●	544 ...	5879-	012	014	016	018				-3
FG ●	534 ...	3879-	012	014	016					-3
FG ●	514 ...	8879-	012	014	016					-3
FG ●	504 ...	4879-	012		016					-3

## 879K

### Long tapered chamfer, torpedo

806 314 299



Shank		L mm									
		10,0	10,0	10,0	10,0	10,0	Shank				
		Order N°						Shank			
FG ●	524 ...	879K-	012	014	016	018	021			-3	
FG ●	544 ...	5879K-	012	014	016	018	021			-3	
FG ●	534 ...	3879K-	012	014	016	018	021			-3	
FG ●	514 ...	8879K-	012	014	016	018	021			-3	
FG ●	504 ...	4879K-	012	014	016					-3	

Match C284K






## 880

### Short parallel chamfer, round

806 314 140




Shank		L mm						
		6,0	6,0	6,0	6,0	6,0	Shank	
		Order N°						
FG●	524 ...	880-	010	012	014	016	018	-3
FG●	544 ...	5880-		012	014	016	018	-3
FG●	534 ...	3880-	010	012	014	016	018	-3
FG●	514 ...	8880-	010	012	014	016	018	-3
FG●	504 ...	4880-	010	012	014	016		-3

## 881

### Parallel chamfer, round

806 314 141




Shank		L mm					
		8,0	8,0	8,0	8,0	Shank	
		Order N°					
FG●	524 ...	881-	010	012	014	016	-3
FG●	544 ...	5881-			014	016	-3
FG●	534 ...	3881-		012	014	016	-3
FG●	514 ...	8881-	010	012	014		-3
FG●	504 ...	4881-		012	014		-3

## 883

### Pointed

806 314 539




Shank		L mm			
		3,0	Shank		
		Order N°			
FG●	524 ...	883	010	-3	
FG●	534 ...	3883	010	-3	

## 885

### Beveled cylinder

806 314 130




Shank		L mm					
		8,0	8,0	8,0	8,0	Shank	
		Order N°					
FG●	524 ...	885-	012		016		-3
FG●	544 ...	5885-		014		018	-3
FG●	534 ...	3885-	012		016		-3
FG●	514 ...	8885-	012		016		-3



# 886

## Long beveled cylinder

806 314 131



Shank

L mm 10,0 10,0



Order N° Shank


FG●	524 ...	886-	014	016	-3
FG●	544 ...	5886-		016	-3



# 888

## Flame, slim version

806 314 496



Shank

L mm 8,0 8,0 8,0 8,0 8,0



Order N° Shank


FG●	524 ...	888-	010	012	014	016	018	-3
FG	544 ...	5888-		012	014	016	018	-3
FG●	534 ...	3888-	010	012	014	016	018	-3
FG●	514 ...	8888-	010	012	014	016	018	-3
FG●	504 ...	4888-	010	012	014	016		-3



# 889

## Needle

806 314 540



Shank

L mm 3,5 4,0 4,0



Order N° Shank


FG●	524 ...	889-	009	010		-3
FG●	544 ...	5889-		010	012	-3
FG●	534 ...	3889-	009	010		-3
FG●	514 ...	8889-	009	010		-3
FG●	504 ...	4889-	009	010		-3



# 899

## Occlusal/palatal reduction

806 314 033



Shank

L mm 6,5 7,0 7,0



Order N° Shank


FG●	524 ...	899-	021	027	031	-3
FG●	544 ...	5899-	021	027	031	-3
FG●	534 ...	3899-	021	027	031	-3



## 909

### Occlusal reduction

806 314 068




Shank		L mm				
		1,1	2,1	2,2		
		Order N°			Shank	
FG ●	524 ...	909-	040	050	055	-3
FG ●	544 ...	5909-	040	050	055	-3
FG ●	534 ...	3909-	040	050	055	-3



## 970

### Occlusal shaping/finishing

806 314 463




Shank		L mm				
		1,6	1,6			
		Order N°		Shank		
FG ●	524 ...	970-	023	027	-3	
FG ●	514 ...	8970-	023	027	-3	



## 998

### Tapered chamfer with guide-pin

806 314 508



Shank		L mm				
		9,0	9,0	9,0		
		Order N°			Shank	
FG ●	544 ...	5998-	016	021	023	-3
FG ●	514 ...	8998-	016	021	023	-3



# Diamond coated burs

## GOLD



### G 3368

Pointed bud

V806 314 257



Shank

L mm 5,0



Order N° Shank

FG●

534 ...

G3368-

023

-3



### G 3847

Tapered shoulder

V806 314 172



Shank

L mm 8,0 8,0



Order N° Shank

FG●

534 ...

G3847-

014

016

-3



### G 3379

Egg/Football

V806 314 277



Shank

L mm 4,2



Order N° Shank

FG●

534 ...

G3379-

023

-3



### G 3811

Barrel

V806 314 038



Shank

L mm 4,2



Order N° Shank

FG●

534 ...

G3811-

033


-3



### G 3830L

Pear long

V806 314 239



Shank

L mm 4,0 5,0 5,0



Order N° Shank

FG●

534 ...

G3830L-

012

014

016


-3



## G 3850

Tapered round, long

V806 314 199




Shank		L mm					10,0	10,0	10,0
		Order N°					Shank		
FG ●	534 ...	G3850-	012	014	016	-3			



## G 3856

Tapered chamfer, round

V806 314 198




Shank		L mm				8,0	8,0	8,0	8,0
		Order N°				Shank			
FG ●	534 ...	G3856-	014	016	018	021	-3		



## G 3878K

Tapered chamfer, torpedo

V806 314 298



Shank		L mm					8,0	8,0	8,0	8,0	8,0
		Order N°					Shank				
FG ●	534 ...	G3878K-	014	016	018	020	022	-3			



## G 3879

Long parallel chamfer, torpedo

V806 314 290




Shank		L mm		10,0
		Order N°		Shank
FG ●	534 ...	G3879-	014	-3



## G 3879K

Long tapered chamfer, torpedo

V806 314 299



Shank		L mm				10,0	10,0	10,0	10,0
		Order N°				Shank			
FG ●	534 ...	G3879K-	012	014	016	018	-3		


# Diamond coated burs

## MICROPREPARATION



### MP 801

Micro round




Shank		L mm			
		Order N°			Shank
FG ●	524 ...	MP801	009	012	-3
FG ●	514 ...	MP8801	009	012	-3



### MP 802

Micro round with collar




Shank		L mm				
		2,5	3,0	3,0		
		Order N°				Shank
FG ●	524 ...	MP802	008	012	014	-3
FG ●	514 ...	MP8802	008	012	014	-3



### MP 830

Micro pear




Shank		L mm				
		2,7	2,7			
		Order N°				Shank
FG ●	524 ...	MP830	009	012	-3	
FG ●	514 ...	MP8830	009	012	-3	



### MP 830RM

Micro pear slim



Shank		L mm				
		2,5	2,5			
		Order N°				Shank
FG ●	524 ...	MP830RM	008	010	-3	
FG ●	514 ...	MP8830RM	008	010	-3	





## MP 835KR

Micro cylinder, rounded edge




Shank		L mm		2,0	
		Order N°		Shank	
FG ●	524 ...	MP835KR	008	-3	
FG ●	514 ...	MP8835KR	008	-3	



## MP 838

Micro cylinder, rounded edge




Shank		L mm		2,6	
		Order N°		Shank	
FG ●	524 ...	MP838	007	-3	
FG ●	514 ...	MP8838	007	-3	



## MP 889

Micro lance



Shank		L mm		2,7 3,9	
		Order N°		Shank	
FG ●	524 ...	MP889	007	009	-3
FG ●	514 ...	MP8889	007	009	-3





# Diamond coated burs

## MINI



### 801M

Round

806 313 001



Shank

L mm

		Order N°		Shank	
FG ●	524 ...	801M-	008	-3	



### 845M

Flat end taper

806 313 168



Shank

L mm

		Order N°		Shank	
FG ●	524 ...	845M-	008	-3	



### 805M

Inverted cone

806 313 010



Shank

L mm

		Order N°		Shank	
FG ●	524 ...	805M-	008	-3	



### 852M

Short pointed

806 313 164



Shank

L mm

		Order N°		Shank	
FG ●	524 ...	852M-	008	-3	



### 835M

Cylinder flat end

806 313 107



Shank

L mm

		Order N°		Shank	
FG ●	524 ...	835M-	008	-3	



# Diamond coated burs

## RAPID GRINDER



### 837S

Rapid grinder, flat end

806 314 A111



Shank		L mm 7,0		
Order N°		Shank		
FG ●	524 ...	837S	014	-3
FG ●	534 ...	3837S	014	-3



### 848A

Rapid grinder, round end

806 314 516



Shank		L mm 10,0		
Order N°		Shank		
FG ●	524 ...	848A	018	-3
FG ●	534 ...	3848A	018	-3



### 842A

Rapid grinder, round end

806 314 510



Shank		L mm 10,0		
Order N°		Shank		
FG ●	524 ...	842A	018	-3
FG ●	534 ...	3842A	018	-3



### 848S

Rapid grinder

806 314 A173



Shank		L mm 10,0 10,0			
Order N°		Shank			
FG ●	524 ...	848S	018	023	-3



### 842S

Rapid grinder, flat end

806 314 A113



Shank		L mm 10,0		
Order N°		Shank		
FG ●	524 ...	842S	018	-3
FG ●	534 ...	3842S	018	-3



### 847S

Rapid grinder

806 314 A172



Shank		L mm 7,0		
Order N°		Shank		
FG ●	524 ...	847S	014	-3




# Diamond coated burs


## TURBO DIAMONDS

### 811T

Occlusal/palatal reduction

T806 314 038




Shank	L mm		4,2		
	Order N°		Shank		
FG ●	534 ...	3811T	018	-3	

### 878KT

Tapered chamfer, torpedo

T806 314 298




Shank	L mm		8,0		
	Order N°		Shank		
FG ●	534 ...	3878KT	014	-3	

### 850T

Tapered round, long

T806 314 199




Shank	L mm		10,0 10,0		
	Order N°		Shank		
FG ●	544 ...	5850T	014	016	-3
FG ●	534 ...	3850T	014	016	-3

### 881T

Parallel chamfer, round

T806 314 141




Shank	L mm		8,0 8,0		
	Order N°		Shank		
FG ●	534 ...	3881T	012	014	-3

### 859T

Needle

T806 314 166




Shank	L mm		10,0		
	Order N°		Shank		
FG ●	544 ...	5859T	018	-3	

### 863T

Long flame

T806 314 250



Shank	L mm				
	Order N°		Shank		
FG ●	534 ...	3863T	014	-3	



# Diamond coated burs

## ZIRCONA



### Z 801

Round

806 314 001



Shank

L mm

Order N°

Shank

FG ●●

535 ...

Z801-

010

-3



### Z 830

Pear

806 314 238



Shank

L mm

Order N°

Shank

FG ●●

535 ...

Z830-

014

-3



### Z 856

Short tapered chamfer, round

806 314 197



Shank

L mm

Order N°

Shank

FG ●●

535 ...

Z856-

014


-3



### Z 880

Short parallel chamfer, round

806 314 140



Shank

L mm

Order N°

Shank

FG ●●

535 ...

Z880-

014

-3





## Z 3379

Egg/Football

806 314 277

Shank

L mm

Order N° Shank

FG ●● 535 ... Z3379- 023 -3



## Z 3801L

Round

806 314 001

Shank

L mm

Order N° Shank

FG ●● 535 ... Z3801L- 014 -3



## Z 3838L

Short cylinder rounded edge

806 314 137

Shank

L mm

Order N° Shank

FG ●● 535 ... Z3838L- 012 -3



## Z 3850

Tapered round, long


806 314 197

Shank

L mm

Order N° Shank

FG ●● 535 ... Z3850- 025 -3



## Z 3860

Short flame

806 314 245

Shank

L mm

Order N° Shank

FG ●● 535 ... Z3860- 014 -3



# Carbide instruments



Carbide	p. 40
Finishing Instrument	p. 48
Crown Cutters	p. 63
Amalgam Remover	p. 65
Intraoral trimming on titanium	p. 67
Surgery	p. 69
Periodontic	p. 75

# Carbide instruments


## CARBIDES



### C1SXM

Round

C1SXM.204. ...



Shank



Order N° Shank

RA ●	C1SXM-	010	012	014	016	018	021	023	027	-2
------	--------	-----	-----	-----	-----	-----	-----	-----	-----	----


500 - 3000 rpm



### C1SXML

Round, long

C1SXM.205. ...



Shank



Order N° Shank

RA L ●	C1SXML-	010	014	018	023	-2
--------	---------	-----	-----	-----	-----	----


500 - 3000 rpm



### C1S

Round

500 314 001 003



Shank



Order N° Shank

FG ●	C1S-	010	012	014	016	018	021	023		-3
------	------	-----	-----	-----	-----	-----	-----	-----	--	----

Max 300 000 rpm

500 204 001 003

RA ●	C1S-	010	012	014	016	018	021	023	027	-2
------	------	-----	-----	-----	-----	-----	-----	-----	-----	----

Max 50 000 rpm






# C1SL

Round, long

500 205 001 003



Shank



Order N°		Shank					
RA L	C1SL-	010	014	018	023	027	-2


Max 50 000 rpm



# C1SNL

Round, extra-long

500 206 001 003



Shank



Order N°		Shank			
RA XL	C1SNL-	006	008	010	-2


1500 - 2000 rpm



# C1

Round

500 314 001 001



Shank



Order N°													Shank	
FG	C1-	005	006	007	008	010	012	014	016	018	021	023	027	-3

Max 300 000 rpm

500 204 001 001

Order N°		Shank										
RA	C1-	007	008	010	012	014	016	018	021	023	027	-2


Max 50 000 rpm



# C1L

Round, long

500 316 001 001



Shank

Order N°								Shank
FGXL	C1L-	010	012	014	016	018	023	-3


Max 300 000 rpm



## C1X

Round, cross

500 314 001 190



Shank



Order N°		Shank						
FG		C1X-	010	012	014	016	018	-3


180 000 - 200 000 rpm



## C1XL

Round, long cross

500 206 001 294



Shank


Order N°		Shank					
RAXL		C1XL-	014	018	023	027	-2



## C2

Inverted cone

500 314 010 006



Shank

L mm		Order N°								Shank	
		0,5	0,9	1,0	1,1	1,2	1,4	1,6	1,7		
FG		C2-	006	007	008	010	012	014	016	018	-3

Max 300 000 rpm

500 204 010 006

RA		C2-				010	012	014	016	018	-2
----	--	-----	--	--	--	-----	-----	-----	-----	-----	----

Max 50 000 rpm






## C7

### Pear

500 314 232 001



#### Shank

L mm	1,7	1,8	1,8	1,9
------	-----	-----	-----	-----



Order N°					Shank
----------	--	--	--	--	-------

FG		C7-	008	010	012	014	-3
----	--	-----	-----	-----	-----	-----	----

Max 300 000 rpm

500 204 232 001

RA		C7-	008	010	012		-2
----	--	-----	-----	-----	-----	--	----


Max 50 000 rpm



## C21

### Cylinder

500 314 107 006



#### Shank

L mm	3,5	4,0	4,0	4,5	4,5	5,0
------	-----	-----	-----	-----	-----	-----



Order N°						Shank
----------	--	--	--	--	--	-------

FG		C21-	008	009	010	012	014	016	-3
----	--	------	-----	-----	-----	-----	-----	-----	----


Max 300 000 rpm



## C21L

### Cylinder, long

500 314 110 006



#### Shank

L mm	5,5	5,5	6,0	6,0
------	-----	-----	-----	-----



Order N°					Shank
----------	--	--	--	--	-------

FG		C21L-	009	010	012	014	-3
----	--	-------	-----	-----	-----	-----	----


Max 300 000 rpm



## C21R

### Cylinder, rond

500 314 137 006



#### Shank

L mm	3,5	4,0	4,5	4,5
------	-----	-----	-----	-----



Order N°					Shank
----------	--	--	--	--	-------

FG		C21R-	008	010	012	014	-3
----	--	-------	-----	-----	-----	-----	----

Max 300 000 rpm

500 204 137 006

RA		C21R-	008	010			-2
----	--	-------	-----	-----	--	--	----


Max 50 000 rpm



## C21RX

Cylinder, rond cross

500 314 140 019



Shank

L mm 5,6 5,6



Order N° Shank

FG		C21RX-	010	012	-3
----	--	--------	-----	-----	----


180 000 - 200 000 rpm



## C23

Tapered

500 314 168 006



Shank

L mm 3,5 4,0 4,0 4,5 4,5



Order N° Shank

FG		C23-	008	009	010	012	014	-3
----	--	------	-----	-----	-----	-----	-----	----

Max 300 000 rpm

500 204 168 006

RA		C23-	008		010	012	014	-2
----	--	------	-----	--	-----	-----	-----	----


Max 50 000 rpm



## C23L

Tapered, long

500 314 171 006



Shank

L mm 5,5 6,0



Order N° Shank

FG		C23L-	010	012	-3
----	--	-------	-----	-----	----


Max 300 000 rpm



## C23R

Tapered, rond

500 314 194 006



Shank

L mm 4,0 4,5 4,5 4,5



Order N°

Shank

FG	C23R-	010	012	014	016	-3
----	-------	-----	-----	-----	-----	----

Max 300 000 rpm

500 204 194 006


RA	C23R-	010	012	014	016	-2
----	-------	-----	-----	-----	-----	----

Max 50 000 rpm

## C23RX

Tapered, rond cross

500 314 197 019



Shank

L mm 5,6 5,6



Order N°

Shank


FG	C23RX-	010	012	-3
----	--------	-----	-----	----

180 000 - 200 000 rpm

## C23SR

Tapered, small rond

500 314 194 019



Shank

L mm 4,0 4,0



Order N°

Shank


FG	C23SR-	008	010	-3
----	--------	-----	-----	----

Max 300 000 rpm

## C25R

Pear, long

500 314 238 006



Shank

L mm 4,0 4,5 4,5



Order N°

Shank


FG	C25R-	010	012	014	-3
----	-------	-----	-----	-----	----

Max 300 000 rpm

## C31

### Cylinder with cross cut

500 314 107 007



Shank

L mm 3,5 4,0 4,5 4,5 5,0



Order N° Shank

FG	C31-	008	010	012	014	016	-3
----	------	-----	-----	-----	-----	-----	----

Max 300 000 rpm

500 204 107 007


RA	C31-	008	010	012	014		-2
----	------	-----	-----	-----	-----	--	----

Max 50 000 rpm

## C31L

### Cylinder, long with cross cut

500 314 110 007



Shank

L mm 5,5 6,0



Order N° Shank


FG	C31L-	010	012	-3
----	-------	-----	-----	----

Max 300 000 rpm

## C31R

### Cylinder, rond with cross cut

500 314 137 007



Shank

L mm 3,5 4,0 4,5 4,5



Order N° Shank


FG	C31R-	008	010	012	014	-3
----	-------	-----	-----	-----	-----	----

Max 300 000 rpm

## C33

Tapered with cross cut

500 314 168 007



Shank

L mm 4,0 4,0 4,5 4,5 5,0



Order N°

Shank

FG		C33-	009	010	012	014	016	-3
----	--	------	-----	-----	-----	-----	-----	----

Max 300 000 rpm

500 204 168 007


RA		C33-			012			-2
----	--	------	--	--	-----	--	--	----

Max 50 000 rpm

## C33L

Tapered, long with cross cut

500 314 171 007



Shank

L mm 5,5 5,5 6,0 6,5



Order N°

Shank


FG		C33L-	009	010	012	016	-3
----	--	-------	-----	-----	-----	-----	----

Max 300 000 rpm

## C33R

Tapered, rond with cross cut

500 314 194 007



Shank

L mm 4,0 4,5



Order N°

Shank

FG		C33R-	010	014	-3
----	--	-------	-----	-----	----

Max 300 000 rpm

# Carbide instruments

## FINISHING INSTRUMENTS



### C41

Round

500 314 001 071



Shank

L mm 1,2 1,6



Order N° Shank

FG ●

C41- 014 018 -3

Max 300 000 rpm up to diam 018

Max 140 000 rpm from diam 023

500 204 001 071

RA ●

C41- 014 018 -2

Max 100 000 rpm up to diam 018

Max 50 000 rpm from diam 023

12 blades 014

18 blades 018

20 blades 023



### C41F

Round

500 314 001 041



Shank

L mm 1,6



Order N° Shank

FG ●

C41F- 018 -3

Max 300 000 rpm

500 204 001 041

RA ●

C41F- 018 -2

Max 50 000 rpm

22 blades





## C44E

Egg/Football

500 314 277 072

Shank




FG ●

L mm 2,9 3,5 4,2

Order N° Shank

C44E- 014 018 023 -3



Max 300 000 rpm

500 204 277 072

RA ●

C44E- 014 018 023 -2

Max 50 000 rpm

12 blades



## C44EF

Egg/Football

500 314 277 042

Shank



FG ●

L mm 3,5

Order N° Shank

C44EF- 018 -3



Max 300 000 rpm

500 204 277 042

RA ●

C44EF- 018 -2

Max 50 000 rpm

16-22 blades



## C44UF

### Egg/Football

500 314 277 032



Shank

L mm 3,5



Order N° Shank

FG

C44UF- 018 -3

Max 300 000 rpm

500 204 277 032

RA

C44UF- 018 -2

Max 50 000 rpm

30 blades



## C45

### Grenade

500 314 274 072



Shank

L mm 4,0



Order N° Shank

FG

C45- 018 -3

Max 300 000 rpm

500 204 274 072

RA

C45- 018 -2

Max 50 000 rpm

12 blades





## C45F

### Grenade

500 314 274 042



#### Shank

L mm 4,0



Order N° Shank

FG ●

C45F- 018 -3

Max 300 000 rpm

500 204 274 042

RA ●

C45F- 018 -2

Max 50 000 rpm


16-22 blades



## C46E

### Flame

500 314 254 072



#### Shank

L mm 2,4 2,9 3,3 3,3 3,6 4,1



Order N° Shank

FG ●

C46E- 009 012 014 018 023 -3

Max 300 000 rpm

500 204 254 072

RA ●

C46E- 010 012 014 -2

Max 50 000 rpm


12 blades



## C46EL

### Flame

500 314 496 071



#### Shank

L mm 5,4 5,4



Order N° Shank

FG ●

C46E- 010 012 -3


Max 300 000 rpm



## C48

### Flame

500 314 243 072



### Shank



L mm	3,4	3,5	3,5	4,0
------	-----	-----	-----	-----

Order N°	Shank				
----------	-------	--	--	--	--

FG ●		C48-	009	010	012	014	-3
------	--	------	-----	-----	-----	-----	----

Max 300 000 rpm

12 blades



## C48F

### Flame

500 314 243 042



### Shank



L mm	3,5
------	-----

Order N°	Shank
----------	-------

FG ●		C48F-	012	-3
------	--	-------	-----	----

Max 300 000 rpm

500 204 243 042

Order N°	Shank	
----------	-------	--

RA ●		C48F-	012	-2
------	--	-------	-----	----

Max 50 000 rpm

16-22 blades



## C48L

### Flame, long

500 314 249 072



### Shank



L mm	8,0
------	-----

Order N°	Shank
----------	-------

FG ●		C48L-	012	-3
------	--	-------	-----	----

Max 300 000 rpm

12 blades

Matches 862 diamond



## C48LF

Flame, long

500 314 249 042



Shank

L mm 8,0



Order N° Shank

FG ●

C48LF- 012 -3

Max 300 000 rpm

16-22 blades

Matches 862 diamond



## C48LUF

Flame, long

500 314 249 032



Shank

L mm 8,0



Order N° Shank

FG ●

C48LUF- 012 -3

Max 300 000 rpm

30 blades

Matches 862 diamond



## C50A

Interproximal, safe end

500 314 474 071



Shank

L mm 8,0



Order N° Shank

FG ●

C50A- 010 -3

Max 300 000 rpm

12 blades



## C50AF

Interproximal, safe end

500 314 474 041



Shank

L mm 8,0



Order N° Shank

FG ●

C50AF- 010 -3

Max 300 000 rpm

16-22 blades





## C332

### Safe end

500 314 699 071

#### Shank



FG ●		C332-	008	-3
------	--	-------	-----	----

Max 300 000 rpm

500 204 699 071

RA ●		C332-	008	-2
------	--	-------	-----	----

Max 50 000 rpm

8 blades



L mm	3,5
------	-----

Order N°	Shank
----------	-------



## C333

### Small point

500 314 159 071

#### Shank



FG ●		C333-	010	-3
------	--	-------	-----	----

Max 300 000 rpm

8 blades



L mm	4,2
------	-----

Order N°	Shank
----------	-------



## C333F

### Small point

500 314 159 041

#### Shank



FG ●		C333F-	010	-3
------	--	--------	-----	----

Max 300 000 rpm

16-22 blades



L mm	4,2
------	-----


Order N°	Shank
----------	-------



## C333UF

### Small point

500 314 159 031



#### Shank



L mm 4,2

Order N° Shank

FG ●

C333UF- 010 -3

Max 300 000 rpm

30 blades



## C334

### Medium point

500 314 164 071



#### Shank



L mm 6,0

Order N° Shank

FG ●

C334- 014 -3

Max 300 000 rpm

8 blades

Matches 8852 diamond



## C334F

### Medium point

500 314 164 041



#### Shank



L mm 6,0

Order N° Shank

FG ●

C334F- 014 -3

Max 300 000 rpm

16-22 blades

Matches 4852 diamond





## C334UF

### Medium point

500 314 164 031



#### Shank



FG		C334UF-	014	-3
----	--	---------	-----	----

L mm	6,0
------	-----

Order N°	Shank
----------	-------

Max 300 000 rpm

30 blades

Matches 852UF diamond



## C335

### Long point

500 314 166 071



#### Shank



FG		C335-	014	-3
----	--	-------	-----	----

L mm	9,0
------	-----

Order N°	Shank
----------	-------

Max 300 000 rpm

8 blades

Matches 8859 diamond



## C335F

### Long point

500 314 166 041



#### Shank



FG		C335F-	014	-3
----	--	--------	-----	----

L mm	9,0
------	-----

Order N°	Shank
----------	-------

Max 300 000 rpm

16-22 blades

Matches 4859 diamond





# C335UF

## Long point

500 314 166 031



### Shank



FG		C335UF-	014	-3
----	--	---------	-----	----

L mm	9,0
------	-----

Order N°	Shank
----------	-------

Max 300 000 rpm

30 blades


Matches 859UF diamond



# C243

## Parallel chamfer, Torpedo

500 314 288 072



### Shank



FG		C243-	010	012	-3
----	--	-------	-----	-----	----

L mm	6,0	6,0
------	-----	-----

Order N°	Shank
----------	-------

Max 300 000 rpm

500 204 288 072

RA		C243-	010	012	-2
----	--	-------	-----	-----	----

Max 100 000 rpm

8 to 10 blades


Matches 877 diamond



# C283

## Parallel chamfer, Torpedo

500 314 289 072



### Shank



FG		C283-	012	014	-3
----	--	-------	-----	-----	----

L mm	8,0	8,0
------	-----	-----

Order N°	Shank
----------	-------

Max 300 000 rpm

500 204 289 072

RA		C283-	012	014	-2
----	--	-------	-----	-----	----

Max 100 000 rpm

8 to 12 blades

Matches 878 diamond



## C282K

### Tapered chamfer, Torpedo

500 314 297 072



Shank

L mm 6,0

2°



Order N° Shank

FG ●

C282K- 014 -3

Max 300 000 rpm

500 204 297 072

RA ●

C282K- 014 -2

Max 100 000 rpm

8 to 10 blades


Matches 877K diamond



## C244K

### Tapered chamfer, Torpedo

500 314 298 072



Shank

L mm 8,0 8,0

2° 2°



Order N° Shank

FG ●

C244K- 014 016 -3

Max 300 000 rpm

500 204 298 072

RA ●

C244K- 014 016 -2

Max 100 000 rpm

10 to 12 blades

Matches 878K diamond






# C283K

Tapered chamfer, Torpedo

500 314 298 072



Shank

L mm 8,0

2°



Order N° Shank

FG ●

C283K- 021 -3

Max 300 000 rpm

500 204 298 072

RA ●

C283K- 021 -2

Max 100 000 rpm

10 to 12 blades

Matches 878K diamond



# C284K

Tapered chamfer, Torpedo

500 314 299 072



Shank

L mm 10,0

2°



Order N° Shank

FG ●

C284K- 018 -3

Max 300 000 rpm

12 blades


Matches 879K diamond



# C244KR

Parallel shoulder, rounded edge

500 314 158 072



Shank

L mm 8,0 8,0



Order N° Shank

FG ●

C244KR- 012 014 -3

Max 300 000 rpm

10 blades

Matches 837KR diamond



## C336

### Tapered shoulder, Torpedo

500 314 198 072



Shank		L mm	8,0		
		Order N°	Shank		
FG ●		C336-	016	-3	

160 000 rpm

12 blades


Matches 847KR diamond



## C375LR

### Tapered chamfer, Round

500 314 198 072



Shank		L mm	8,0	8,0	8,0	8,0	8,0
 			2°	2°	2°	2°	2°
Order N°		Shank					
FG ●		C375LR-	012	014	016	018	023
							-3

Max 300 000 rpm

12 blades

Matches 856 diamond



## C207

### End cutting bur

500 314 150 001



Shank				
		Order N°	Shank	
FG		C207-	012	-3

Matches 839 diamond






# C390

## Grenade

500 314 274 072



### Shank



L mm 3,5 3,5 3,7

Order N° Shank

FG ● C390- 014 016 018 -3

Max 300 000 rpm

500 204 274 072

RA ● C390- 014 016 018 -2

18 000 - 30 000 rpm


12 blades



# C390UF

## Grenade

500 314 274 032



### Shank



L mm 3,5

Order N° Shank

FG ● C390- 016 -3

Max 300 000 rpm

500 204 274 032

RA ● C390- 016 -2

18 000 - 30 000 rpm

30 blades






# C246

## Needle/Flame

500 314 495 071



### Shank



Order N°		Shank						
FG ●		C246-	008	009	010	012	014	-3

Max 300 000 rpm

500 204 495 071

RA ●		C246-				012		-2
------	--	-------	--	--	--	-----	--	----

Max 100 000 rpm


12 blades



# C446

## Square

500 314 467 211



### Shank



Order N°		Shank				
FG		C446-	010	012	014	-3

# C97

## Three-edges wedge

500 314 468 373



### Shank



Order N°		Shank		
FG		C97-	010	-3

# Carbide instruments

## CROWN CUTTERS



### C34


For metal crowns

500 314 138 293

Shank



FG●



L mm 2,0 2,0

Order N° Shank

C34- 010 012 -3

180 000 - 200 000 rpm

Cross-cut



### C36R


For metal crowns

500 314 139 293

Shank



FG●



L mm 4,0 4,5

Order N° Shank

C36R- 010 012 -3

180 000 - 200 000 rpm

Cross-cut



### C21RX


For metal crowns

500 314 140 019

Shank



FG



L mm 5,6 5,6

Order N° Shank

C21RX- 010 012 -3

180 000 - 200 000 rpm


Double-cut



## C23RX

For metal crowns

500 314 197 019



Shank

L mm 5,6 5,6



Order N° Shank

FG		C23RX-	010	012	-3
----	--	--------	-----	-----	----

180 000 - 200 000 rpm



## C21A

For metal crowns

500 314 137 006



Shank

L mm 4,0



Order N° Shank

FG		C21A-	012	-3
----	--	-------	-----	----


Single-use



## C4CC

For metal, titanium and low-fused ceramic crowns

500 314 138 019



Shank

L mm 2,1 2,1



Order N° Shank

FG ●●		C4CC-	010	012	-3
-------	--	-------	-----	-----	----

120 000 - 160 000 rpm



## C4CCL

For metal, titanium and low-fused ceramic crowns

500 314 139 019



Shank

L mm 3,5



Order N° Shank

FG ●●		C4CCL-	012	-3
-------	--	--------	-----	----

120 000 - 160 000 rpm



## C4T

For metal crowns

500 314 139 008



Shank

L mm 4,0



Order N° Shank

FG		C4T-	012	-3
----	--	------	-----	----

60 000 - 160 000 rpm



# Carbide instruments


## AMALGAM REMOVER



### C32HG

Amalgam remover

500 314 137 018 ...



Shank

L mm 4,2



Order N° Shank

FG ●●

C32HG- 012 -3

180 000 - 200 000 rpm


Medium double-cut with overflow cut



### CTC7X

Amalgam & composite remover

500 314 232 019 ...



Shank

L mm 1,8



Order N° Shank

FG

CTC7X- 012 -3

160 000 rpm



### CTC7XL

Amalgam & composite remover

500 314 234 019 ...



Shank

L mm 4,2 4,4



Order N° Shank

FG

CTC7XL- 012 014 -3


160 000 rpm



## CTC21RX

Amalgam & composite remover

500 314 137 019 ...



Shank		L mm		4,2	4,2
		Order N°		Shank	
FG		CTC21RX-	012	014	-3

160 000 rpm



## CTC23RX

Amalgam & composite remover

500 314 194 019 ...



Shank		L mm		4,2
		Order N°		Shank
FG		CTC23RX-	012	-3

160 000 rpm

# Carbide instruments

## INTRAORAL TRIMMING ON TITANIUM



### C275RX

500 314 141 080 ...

Shank



FG

L mm 8,0

Order N° Shank

C275RX- 012 -3



100 000 rpm



### C283RX

500 314 289 080 ...

Shank



FG

L mm 8,0

Order N° Shank

C283RX- 012 -3



100 000 rpm



### C375RX

500 314 198 080 ...

Shank



FG

L mm 8,0

Order N° Shank

C375RX- 016 -3



100 000 rpm



## C379RX

500 314 277 080 ...



Shank

L mm 8,0



Order N° Shank

FG		C379RX-	023	-3
----	--	---------	-----	----

100 000 rpm



## C48RX

500 314 289 080 ...



Shank

L mm 8,0



Order N° Shank

FG		C48RX-	012	-3
----	--	--------	-----	----

100 000 rpm



# Carbide instruments

## SURGERY



### CTC135

Soft tissue trimmer



Shank		L mm	5,0		
		Order N°	Shank		
FG		CTC135-	016	-3	

300 000 - 500 000 rpm



### CTC46

Soft tissue trimmer



Shank		L mm	8,0		
		Order N°	Shank		
FG		CTC46-	016	-3	


300 000 - 500 000 rpm



### C31RL

Bone cutter

500 316 140 107 ...



Shank		L mm	6,0	6,0	
		Order N°	Shank		
FG		C31RL-	010	012	-3

Max 100 000 rpm



## C33RL

### Bone cutter

500 316 140 107 ...



Shank

L mm 6,0 6,0 6,0



Order N° Shank

FG	C33RL-	014	016	018	-3
----	--------	-----	-----	-----	----

Max 100 000 rpm



## C700XXL

### Bone cutter

500 314 415 296 ...



Shank

L mm 5,0



Order N° Shank

FG	C700XXL-	011	-3
----	----------	-----	----

15 000 - 60 000 rpm



## C269

### Bone cutter

500 314 199 295 ...



Shank

L mm 11,0



Order N° Shank

FG	C269-	016	-3
----	-------	-----	----

500 204 199 295 ...

RA	C269-	016	-2
----	-------	-----	----

500 104 199 295 ...

HP	C269-	016	-104
----	-------	-----	------

Max 3000 rpm






## C151

### Bone cutter, Zekrya

500 316 199 336 ...



Shank

L mm 11,0 11,0



Order N° Shank

FG		C151P-	016	-3
----	--	--------	-----	----

500 317 199 336 ...

FG		C151PL-	016	-3
----	--	---------	-----	----

500 204 199 336 ...

RA		C151-	012	-2
----	--	-------	-----	----

500 104 199 336 ...

HP		C151-	012	-1
----	--	-------	-----	----

Max 160 000 rpm



## C152

### Bone cutter, Endo-Z

500 314 210 295 ...



Shank

L mm 9,0



Order N° Shank

FG		C152-	014	-3
----	--	-------	-----	----

500 315 210 295 ...

FG		C152L-	014	-3
----	--	--------	-----	----

Max 120 000 rpm



## C269GK

### Bone cutter, Endo-Z

500 314 219 295 ...



Shank

L mm 9,0



Order N° Shank

FG		C269GK-	016	-3
----	--	---------	-----	----

Max 120 000 rpm





# C161

## Bone cutter, Lindemann

500 314 408 295 ...



Shank		L mm	9,0		
		Order N°	Shank		
FG		C161-	016	-3	

500 104 408 295 ...

HP		C161-	016	-1	
----	--	-------	-----	----	--

Max 3000 rpm



# C162

## Bone cutter, Lindemann

500 314 408 297 ...



Shank		L mm	9,0		
		Order N°	Shank		
FG		C162-	016	-3	

500 204 408 297 ...

RA		C162-	016	-2	
----	--	-------	-----	----	--

500 104 408 297 ...

HP		C162-	016	-1	
----	--	-------	-----	----	--

Max 3000 rpm



# C162A

## Bone cutter, Lindemann

500 314 408 298 ...



Shank		L mm	9,0		
		Order N°	Shank		
FG		C162A-	016	-3	

500 204 408 298 ...

RA		C162A-	016	-2	
----	--	--------	-----	----	--

500 104 408 298 ...

HP		C162A-	016	-1	
----	--	--------	-----	----	--

Max 3000 rpm






# C254

## Bone cutter

500 316 254 297 ...



Shank

L mm 6,0 6,0



Order N° Shank

FG		C254-	012		-316
----	--	-------	-----	--	------

500 317 254 297 ...

FG		C254-	012		-317
----	--	-------	-----	--	------

500 205 254 297 ...

RA		C254-	012	016	-205
----	--	-------	-----	-----	------

500 104 254 297 ...

HP		C254-	012	016	-104
----	--	-------	-----	-----	------

Max 160 000 rpm



# C408S

## Bone cutter

500 316 408S 297 ...



Shank

L mm 10,0



Order N° Shank

FG		C408S-	016		-316
----	--	--------	-----	--	------

500 204 408S 297 ...

RA		C408S-	016		-205
----	--	--------	-----	--	------

500 104 408S 297 ...

HP		C408S-	016		-104
----	--	--------	-----	--	------

500 106 408S 297 ...

HP		C408S-	016		-106
----	--	--------	-----	--	------



# C409S

## Bone cutter

500 205 409S 297 ...



Shank

L mm 10,0



Order N° Shank

RA		C409S-	016	-205
----	--	--------	-----	------

500 104 409S 297 ...

HP		C409S-	016	-104
----	--	--------	-----	------

500 106 409S 297 ...


HP		C409S-	016	-106
----	--	--------	-----	------



# C141F

## Bone cutter round

500 104 001 251



Shank



Order N° Shank

HP		C141F-	010	014	018	023	027	031	035	040	-104
----	--	--------	-----	-----	-----	-----	-----	-----	-----	-----	------

500 205 001 251

RA		C141F-	010	014	018	023	027	031	035	040	-205
----	--	--------	-----	-----	-----	-----	-----	-----	-----	-----	------


Max 3000 rpm



# C141AL

## Bone cutter round

500 104 001 298



Shank



Order N° Shank

HP		C141AL-	023	027	031	035	-104
----	--	---------	-----	-----	-----	-----	------

500 205 001 298

RA		C141AL-	023	027	031	035	-205
----	--	---------	-----	-----	-----	-----	------

Max 3000 rpm

# Carbide instruments

## PERIODONTIC



### C407

500 204 258 381 ...



Shank

L mm 2,9

Order N° Shank

RA		C407-	012	-2
----	--	-------	-----	----

Max 12 000 rpm

Root smoothing



### C407L

500 204 259 381 ...



Shank

L mm 2,9

Order N° Shank

RA		C407L-	012	-2
----	--	--------	-----	----

Max 12 000 rpm

Root smoothing

# Polishing instruments





# Polishing instruments

## SILICONE POLISHERS

### Universal Polishers

Medium



Shank

Shank

RA		Order N°	G1	G5	G5B	G4	G2	G7	G8	G9	-2

7 000 to 10 000 rpm

### Ceramic Polishers

Coarse



Shank

Shank

RA		Order N°	KV2B	KV7	-2

7 000 to 15 000 rpm

### Ceramic Polishers

Medium



Shank

Shank

RA		Order N°	KN2B	KN7	-2

7 000 to 15 000 rpm

### Ceramic Polishers

Fine



Shank

Shank

RA		Order N°	KS2B	KS7	-2

# Polishing instruments

## PROPHYLAXIS POLISHERS

### Prophylaxis



Shank

Shank



RA

Order N°

PR7L

-2

5 000 to 12 000 rpm

### Prophylaxis



Shank

Shank



RA

Order N°

PR7

-2

5 000 to 12 000 rpm

### Prophylaxis



Shank

Shank



RA

Order N°

PR111

-2

5 000 to 12 000 rpm

### Prophylaxis



Shank

Shank



RA

Order N°

PR5

-2

5 000 to 12 000 rpm

## Prophylaxis



Shank		Shank		
RA		Order N°	PR5K	-2

5 000 to 12 000 rpm

## Prophylaxis



Shank		Shank		
RA		Order N°	PR112	-2

5 000 to 12 000 rpm

## Prophylaxis



Shank		Shank		
RA		Order N°	B1-045	-2

5 000 rpm with polishing paste

## Prophylaxis



Shank		Shank		
RA		Order N°	B12-045	-2


5 000 rpm with polishing paste





## Prophylaxis





Shank		Shank		
				
RA		Order N°	B12-065	-2

5 000 rpm with polishing paste

## Prophylaxis





Shank		Shank		
 				
RA		Order N°	DB1-045	-2

6 000 rpm

## Prophylaxis





Shank		Shank		
 				
RA		Order N°	DB12-045	-2


6 000 rpm

## Prophylaxis



Shank		Shank		
 				
RA		Order N°	DB12-065	-2

6 000 rpm



# Polishing instruments

## SYNTHETIC RUBBER POLISHERS

### Gold, Amalgam and composite

Removing and contouring



Shank

Shank



RA		Order N°	N1L	N1S	N1	N2	N2W	N6			-2
FG		Order N°							N1	N2	-3

### Smoothing



Shank

Shank



RA		Order N°	XX1L	XX1S	XX1	XX2	XX2W	XX6			-2
FG		Order N°							XX1	XX2	-3

### High gloss polishing



Shank

Shank



RA		Order N°	X1L	X1S	X1	X2	X2W	X6			-2
FG		Order N°							X1	X2	-3

## Composite

Smoothing and pre-polishing



Shank

Shank



RA		Order N°	COM7	COM2	COM3	COM6	COM18	-2
----	--	----------	------	------	------	------	-------	----

Polishing



Shank

Shank



RA		Order N°	COM7F	COM2F	COM3F	COM6F	COM18F	-2
----	--	----------	-------	-------	-------	-------	--------	----

## Bracket Remover



Shank

Shank



RA		Order N°	BR3	BR7	-2
----	--	----------	-----	-----	----

5 000 to 7 000 rpm

# Polishing instruments

## DIAMOND POLISHERS

### Ceramic and metal Polishers

Removing and shaping



Shank

Shank



RA		Order N°	6310	6362	6360S	6360	6302-055	6318	6309-120	-2
FG		Order N°		6332						-3

7 000 to 12 000 rpm

### Ceramic and metal Polishers

Smoothing



Shank

Shank



RA		Order N°	6210	6262	6260S	6260	6202-055	6218	6209-120	-2
FG		Order N°		6262						-3

7 000 to 12 000 rpm

### Ceramic and metal Polishers

High gloss polishing



Shank

Shank



RA		Order N°	6110	6162	6160S	6160	6102-055	6118	6109-120	-2
FG		Order N°		6162						-3

7 000 to 12 000 rpm

# All composite materials Polishers

Smoothing and pre-polishing



Shank

Shank



RA	Order N°	5110	5162	5160	5102-055	5118	5107	-2
FG	Order N°		5162					-3

3 000 to 8 000 rpm

# All composite materials Polishers

High gloss polishing



Shank

Shank



RA	Order N°	5210	5262	5260	5202-055	5218	5207	-2
FG	Order N°		5262					-3

3 000 to 8 000 rpm

# All composite materials Polishers

One-step



Shank

Shank



RA	Order N°	DIA COM2	DIA COM3	DIA COM6	DIA COM7	-2
----	----------	----------	----------	----------	----------	----

5 000 to 7 000 rpm

# Lithium Disilicate, Feldspathic Porcelain

## Medium



Shank

Shank



RA		Order N°	DIA	DIA	DIA	DIA	-2
			LI2	LI3	LI6	LI18	

7 000 to 10 000 rpm

## Fine



Shank

Shank



RA		Order N°	DIA	DIA	DIA	DIA	-2
			LI2f	LI3f	LI6f	LI18f	

7 000 to 10 000 rpm

# Zirconia, Alumina

## Medium



Shank

Shank



RA		Order N°	DIA	DIA	DIA	DIA	-2
			ZIR2	ZIR3	ZIR6	ZIR18	

7 000 to 10 000 rpm

## Fine



Shank

Shank



RA		Order N°	DIA	DIA	DIA	DIA	-2
			ZIR2f	ZIR3f	ZIR6f	ZIR18f	

7 000 to 10 000 rpm

# Polishing instruments

## CERAMIC ABRASIVES

### Abrasives



Shank

Shank



RA		Order N°	K601-030	K652-027	K661-025	K662-035	-2
FG		Order N°	K601-030	K652-027	K661-025	K662-035	-3

FG 15 000 rpm

RA 7 000 rpm

# Polishing instruments

## COMPOSITE/CERAMIC ABRASIVES

### Arkansas



Shank

Shank



		Shank				
RA		Order N°	ARKI-025	ARKI-031	ARK6-030	-2
FG		Order N°	ARKI-025	ARKI-031	ARK6-030	-3

FG 15 000 rpm

RA 7 000 rpm



# Polishing instruments

## FLEXIBLE POLISHING DISC



### Flexible disc

Pack: 100 pieces

Mandrel: 395-2



Order N°	FD10	8FD10	4FD10	FD10UF
Grit	coarse	medium	fine	super-fine
Sizes (mm)	10 x 0,21	10 x 0,17	10 x 0,14	10 x 0,11

Pack: 100 pieces

Mandrel: 395-2



Order N°	FD14	8FD14	4FD14	FD14UF
Grit	coarse	medium	fine	super-fine
Sizes (mm)	14 x 0,21	14 x 0,17	14 x 0,14	14 x 0,11

Blue = Removing and contouring

Red = Smoothing and pre-polishing

Yellow = Polishing

White = High gloss

7 000 to 10 000 rpm

### Starter kit for flexible disc

Pack: 25 pieces of each grit (100 discs)

One mandrel: 395-2



Order N°	Flexi D10	Flexi D14

7 000t to 10 000 rpm



## Diamond Flexible disc

Pack: 100 pieces

Mandrel: 395-2



Order N°	5FDIA10	3FDIA10	FDIA10	8FDIA10
Grit	coarse	medium	fine	super-fine
Sizes (mm)	10 x 0,21	10 x 0,17	10 x 0,14	10 x 0,11

Pack: 100 pieces

Mandrel: 395-2



Order N°	5FDIA14	3FDIA14	FDIA14	8FDIA14
Grit	coarse	medium	fine	super-fine
Sizes (mm)	14 x 0,21	14 x 0,17	14 x 0,14	14 x 0,11

*Black = Removing and contouring*

*Blue = Smoothing and pre-polishing*

*Pink = Polishing*

*Grey = High gloss*

*7 000 to 10 000 rpm*

## Starter kit for flexible disc

Pack: 25 pieces of each grit (100 discs)

One mandrel: 395-2



Order N°	Flexi DIA10	Flexi DIA14

*7 000 to 10 000 rpm*

# Polishing instruments

## DIAMOND STRIPS

### Diamond Finishing Strips

Pack: 10 pieces

Non-serrated diamond strip



Order N°	DS	8DS	4DS
L (mm)	147	147	147
Grit	medium	fine	super-fine
Thickness (mm)	0,13	0,10	0,08
Width (mm)	2,5	2,5	2,5

Pack: 10 pieces


Non-serrated diamond strip



Order N°	DSL	8DSL	4DSL
L (mm)	147	147	147
Grit	medium	fine	super-fine
Thickness (mm)	0,13	0,10	0,08
Width (mm)	3,75	3,75	3,75


Pack: 10 pieces

Serrated diamond strip



Order N°	8DSX	4DSX
L (mm)	147	147
Grit	fine	super-fine
Thickness (mm)	0,10	0,08
Width (mm)	3,75	3,75

# Instruments sets





# Instruments sets


## Set/Crown preparation

KITPREPACO




## Set/Implant abutments molding

KITRX




# Set/Fast, low vibration excavating with clear view

KITC1SXM




# Set/Amalgam and gold

KITPOLAMAL




# Set/Composite


## KITFINCOMPO



## KITPOLCOMP



## KITDIAGOLD





## Set/Lithium disilicate

### KITDIALI




## Set/Zircone

### KITDIAZIR




## Set/Veneer

### KITFACETTE




# Set/Ceramic

KITDIAPOLCER



# Set/Corrective work on prothese

KITRETPROTHESE



# TABLE OF CONTENTS

Diamond.....	p.06
Gold.....	p.08
Micropreparation.....	p.29
Mini.....	p.31
Rapid Grinder.....	p.33
Turbo Diamonds.....	p.34
Zircona.....	p.35
Carbide.....	p.38
Finishing Instrument.....	p.40
Crown Cutters.....	p.48
Amalgam Remover.....	p.63
Intraoral trimming on titanium.....	p.65
Surgery.....	p.68
Periodontic.....	p.69
Silicone Polishers.....	p.75
Prophylaxis Polishers.....	p.78
Synthetic rubber Polishers.....	p.79
Diamond Polishers.....	p.82
Ceramic Abrasives.....	p.84
Composite/Ceramic Abrasives.....	p.87
Flexible Polishing Discs.....	p.88
Diamond Strips.....	p.89
Crown preparation.....	p.91
Implant abutments molding.....	p.94
Fast, low vibration excavating with clear view.....	p.95
Amalgam and gold.....	p.95
Composite.....	p.96
Lithium disilicate.....	p.97
Zicrone.....	p.97
Veneer.....	p.97
Ceramic.....	p.98
Corrective work on prothese.....	p.98



Rue des Anciens Étangs 42  
1190 Bruxelles

Tel: +32 2 343 65 30  
Fax: +32 2 346 28 72

[info@dumont-instruments.com](mailto:info@dumont-instruments.com)  
[www.dumont-instruments.com](http://www.dumont-instruments.com)

***Dumont***  
**INSTRUMENTS & CO**  
Quality since 1978